

# CERTIFICATE OF ANALYSIS

# NorthBio

FETAL BOVINE SERUM  
Cell Culture Grade  
Canadian Origin

Lot # : SF07012  
Catalog # : NBSF-702  
Lot Expiry : May 2015

Test	Test Method	Specification	Result
<i>Sterility and Endotoxin</i>			
Sterility	USP <71>	No evidence of microbial growth	No evidence of microbial growth
Mycoplasma	Large Volume Direct Culture	Not Detected	Not Detected
Mycoplasma	FDA Points to Consider 1993	Not Detected	Not Detected
Endotoxin	USP <85> LAL Gel Clot	< 10 EU/ml	0.077 EU/ml
<i>Viral Safety Testing</i>			
BVD	9 CFR 113.53	Not Detected	Not Detected
HA	9 CFR 113.53	Not Detected	Not Detected
CPE	9 CFR 113.53	Not Detected	Not Detected
Bovine Parvo	9 CFR 113.53	Not Detected	Not Detected
Bovine Adeno 3	9 CFR 113.53	Not Detected	Not Detected
Bovine Adeno 5	9 CFR 113.53	Not Detected	Not Detected
BRSV	9 CFR 113.53	Not Detected	Not Detected
REO Virus	9 CFR 113.53	Not Detected	Not Detected
Rabies	9 CFR 113.53	Not Detected	Not Detected
Blue Tongue	9 CFR 113.53	Not Detected	Not Detected
BVD Antibody	Serum Neutralization	As Reported	1 : 3
<i>Biochemistry</i>			
Total Protein	Biuret	30 – 45 g/L	39 g/L
Albumin	Gel Electrophoresis	As Reported	29.5 g/L
Globulins	Gel Electrophoresis	As Reported	9.5 g/L
Bilirubin	Diazo	As Reported	2.9 µmol/L
ALT	Int'l Fed. Clinical Chemistry	As Reported	6 U/L
AST	Int'l Fed. Clinical Chemistry	As Reported	28 U/L
LDH	Gay's Lactate	As Reported	391 mmol/L
Cholesterol	Cholesterol esterase/oxidase/perox	As Reported	0.4 mmol/L
Glucose	Glucose oxidase/peroxidase	As Reported	2.1 mmol/L
Urea	Urease/glutamate dehydrogenase	As Reported	4.3 mmol/L
Creatinine	Alkaline Picric Acid	As Reported	160 µmol/L
Sodium	Ion Selective Electrode	As Reported	129 mmol/L
Potassium	Ion Selective Electrode	As Reported	18.4 mmol/L
Calcium	Cresolphthalein Complexone	As Reported	2.3 mmol/L
Phosphorus	Ammonium Molybdate	As Reported	2.7 mmol/L
Osmolality	Freezing Point Depression	280 – 340 mOsm/Kg	322 mOsm/Kg
Serum Identity	Ouchterlony	Bovine	Bovine
IgG	Elisa	< 250 mg/mL	86 mg/L
Hemoglobin	Spectrophotometric	≤ 20.0 mg/dL	5.1 mg/dL
pH	Potentiometric	6.80 – 8.20	7.07
Specific Gravity	Weight / Volume	≥ 1.010 g/mL	1.021 g/mL
ALP	Colorometric	As Reported	139 U/L
GGT	Enzymatic Colorometric	As Reported	3 U/L
Chloride	Indirect ISE	As Reported	91 mmol/L
Uric Acid	Enzymatic Colorometric	As Reported	62 µmol/L
Triglyceride	Enzymatic	As Reported	0.38 mmol/L
<i>Cell Culture Performance</i>			
Plating Efficiency	Diploid Cell Line	Pass	Pass
Population Dblg	Diploid Cell Line	Pass	Pass
Cell Toxicity	Diploid Cell Line	None	None
Growth Promotion	Diploid Cell Line	≥ 75 % of in house control	96.0 %
Appearance	Visual	Clear yellow-amber	Complies

### Statement of Origin

THIS SERUM HAS BEEN DERIVED FROM BLOOD ASEPTICALLY COLLECTED IN APPROVED ABATOIRS. THE SERUM WAS ASEPTICALLY PROCESSED AND FILTER STERILIZED IN CANADA IN A FACILITY REGISTERED WITH MEDICAL DEVICES (U.S.) FOLLOWING GUIDELINES OF F.D.A. AND H.P.F.B. (CANADA). THE COUNTRY OF ORIGIN OF THE CRUDE BLOOD (DONOR COW ) IS CANADA.

JULY 28, 2010

DATE

C. BROHMAN, QUALITY ASSURANCE

NORTHBIO INC.  
TORONTO ON CANADA  
M5R 2E3

TEL 1 888 428-9775  
FAX 1 888 857-9775  
EM: [INFO@NORTHBIO.COM](mailto:INFO@NORTHBIO.COM)

